

NRS (Asparaginyl-tRNA synthetase)

Polyclonal Antibody

Cat.NO.: PA01366

3th Edition

Description:Asparaginyl-tRNA synthetase catalyzes the ligation of asparagine with asparaginyl-tRNA using ATP.

Antigen:Full length native recombinant human NRS.

Form:Anti-NRS rabbit serum containing 0.05 % sodium azide for preservative.

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:This antibody reacts specifically with 63kDa asparaginyl tRNA synthetase protein.