

Anti-Human/Mouse/Rat Acetyl-CREBBP (Lys1535) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01416

3th Edition

Description:Acetylates histones, giving a specific tag for transcriptional activation. Also acetylates non-histone proteins, like NCOA3 coactivator. Binds specifically to phosphorylated CREB and enhances its transcriptional activity toward cAMP-responsive genes. Acts as a coactivator of ALX1 in the presence of EP300.

Antigen:Synthesized peptide derived from human CBP around the acetylation site of K1535

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: