

Anti-Human/Mouse/Rat Phospho-AKT2 (Ser474) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01499

3th Edition

Description: General protein kinase capable of phosphorylating several known proteins.

Antigen: Synthesized peptide derived from human Akt2 around the phosphorylation site of Ser474

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in all human cell types so far analyzed.