

Anti-Human/Mouse/Rat Phospho-JUN/JUND (Ser73/100) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01501

3th Edition

Description: Transcription factor that recognizes and binds to the enhancer heptamer motif 5'-TGA[CG]TCA-3'. Promotes activity of NR5A1 when phosphorylated by HIPK3 leading to increased steroidogenic gene expression upon cAMP signaling pathway stimulation.

Antigen: Synthesized peptide derived from human AP-1 around the phosphorylation site of Ser73/100

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: