

**Anti-Human/Mouse/Rat Phospho-MAP3K5 (Ser966) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01505**

---

3th Edition

**Description:**Component of a protein kinase signal transduction cascade. Phosphorylates and activates MAP2K4 and MAP2K6, which in turn activate the JNK and p38 MAP kinases, respectively. Overexpression induces apoptotic cell death.

**Antigen:**Synthesized peptide derived from human ASK 1 around the phosphorylation site of Ser966

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Abundantly expressed in heart and pancreas.