

Anti-Human/Mouse/Rat Phospho-MAP3K5 (Ser966) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01505

3th Edition

Description:Component of a protein kinase signal transduction cascade. Phosphorylates and activates MAP2K4 and MAP2K6, which in turn activate the JNK and p38 MAP kinases, respectively. Overexpression induces apoptotic cell death.

Antigen:Synthesized peptide derived from human ASK 1 around the phosphorylation site of Ser966

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Abundantly expressed in heart and pancreas.