

Anti-Human/Mouse/Rat Phospho-ATF2 (Thr71) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01508

3th Edition

Description: Transcriptional activator, probably constitutive, which binds to the cAMP-responsive element (CRE) (consensus: 5'-GTGACGT[AC][AG]-3'), a sequence present in many viral and cellular promoters. Interaction with JUN redirects JUN to bind to CREs preferentially over the 12-O-tetradecanoylphorbol-13-acetate response elements (TRES) as part of an ATF2/JUN complex.

Antigen: Synthesized peptide derived from human ATF-2 around the phosphorylation site of Thr71

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Abundant expression seen in the brain.