

Anti-Human/Mouse/Rat Phospho-BAD (Ser136) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01511

3th Edition

Description: Promotes cell death. Successfully competes for the binding to Bcl-X(L), Bcl-2 and Bcl-W, thereby affecting the level of heterodimerization of these proteins with BAX. Can reverse the death repressor activity of Bcl-X(L), but not that of Bcl-2 (By similarity). Appears to act as a link between growth factor receptor signaling and the apoptotic pathways.

Antigen: Synthesized peptide derived from human Bad around the phosphorylation site of Ser136

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in a wide variety of tissues.