

Anti-Human/Mouse Phospho-BCR (Tyr177) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01515

3th Edition

Description: GTPase-activating protein for RAC1 and CDC42. Promotes the exchange of RAC or CDC42-bound GDP by GTP, thereby activating them. Displays serine/threonine kinase activity.

Antigen: Synthesized peptide derived from human Bcr around the phosphorylation site of Tyr177

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: