

**Anti-Human/Mouse Phospho-BCR (Tyr177) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01515**

---

3th Edition

**Description:** GTPase-activating protein for RAC1 and CDC42. Promotes the exchange of RAC or CDC42-bound GDP by GTP, thereby activating them. Displays serine/threonine kinase activity.

**Antigen:** Synthesized peptide derived from human Bcr around the phosphorylation site of Tyr177

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**