

**Anti-Human/Mouse/Rat/Monkey Phospho-CDK1/2/3 (Thr14) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01524**

---

3th Edition

**Description:** Plays a key role in the control of the eukaryotic cell cycle. It is required in higher cells for entry into S-phase and mitosis. p34 is a component of the kinase complex that phosphorylates the repetitive C-terminus of RNA polymerase II.

**Antigen:** Synthesized peptide derived from human Cdk1/2/3 around the phosphorylation site of Thr14

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** breast cancer tissues, cancer cell lines and glioblastoma tissue.