

Anti-Human/Mouse/Rat Phospho-CFL1 (Ser3) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01532

3th Edition

Description: Controls reversibly actin polymerization and depolymerization in a pH-sensitive manner. It has the ability to bind G- and F-actin in a 1:1 ratio of cofilin to actin. It is the major component of intranuclear and cytoplasmic actin rods.

Antigen: Synthesized peptide derived from human Cofilin Ser3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely distributed in various tissues.