

Anti-Human/Mouse/Rat Phospho-Cyclin B1 (Ser126) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01538

3th Edition

Description:Essential for the control of the cell cycle at the G2/M (mitosis) transition.

Antigen:Synthesized peptide derived from human Cyclin B1 around the phosphorylation site of Ser126

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: