

**Anti-Human/Mouse/Rat Phospho-FOXO3 (Ser253) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01559**

---

3th Edition

**Description:** Transcriptional activator which triggers apoptosis in the absence of survival factors, including neuronal cell death upon oxidative stress. Recognizes and binds to the DNA sequence 5'-[AG]TAAA[TC]A-3'.

**Antigen:** Synthesized peptide derived from human FoxO3A around the phosphorylation site of Ser253

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitous.