

Anti-Human/Mouse/Rat Phospho-FOXO3 (Ser253) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01559

3th Edition

Description: Transcriptional activator which triggers apoptosis in the absence of survival factors, including neuronal cell death upon oxidative stress. Recognizes and binds to the DNA sequence 5'-[AG]TAAA[TC]A-3'.

Antigen: Synthesized peptide derived from human FoxO3A around the phosphorylation site of Ser253

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitous.