

Anti-Human/Mouse/Rat Phospho-JUNB (Ser259) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01585

3th Edition

Description: Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'.

Antigen: Synthesized peptide derived from human Jun B around the phosphorylation site of Ser259

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: