

Instruction manual FOR RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

Anti-Human/Mouse/Rat/Monkey Phospho-RAF1 (Ser259) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01625

3th Edition

Description:Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. Part of the Ras-dependent signaling pathway from receptors to the nucleus. Protects cells from apoptosis mediated by STK3.

Antigen: Synthesized peptide derived from human Raf-1 around the phosphorylation site of Ser259

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: In skeletal muscle, isoform 1 is more abundant than isoform 2.

1/1