

Anti-Human/Mouse/Rat Phospho-STAT3 (Tyr705) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01632

3th Edition

Description:None

Antigen:Synthesized peptide derived from human Stat3 around the phosphorylation site of Tyr705

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas.