

Anti-Human/Mouse Phospho-EIF4B (Ser422) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01646

3th Edition

Description: EIF4B (eukaryotic translation initiation factor 4B) is a protein-coding gene. GO annotations related to this gene include RNA binding and translation initiation factor activity. An important paralog of this gene is EIF4H.

Antigen: Synthesized peptide derived from human eIF4B around the phosphorylation site of Ser422

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: