

**Anti-Human/Mouse Phospho-EIF4B (Ser422) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01646**

---

3th Edition

**Description:** EIF4B (eukaryotic translation initiation factor 4B) is a protein-coding gene. GO annotations related to this gene include RNA binding and translation initiation factor activity. An important paralog of this gene is EIF4H.

**Antigen:** Synthesized peptide derived from human eIF4B around the phosphorylation site of Ser422

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**