

Anti-Human/Mouse/Rat Phospho-Cyclin D1 (Thr286) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01657

3th Edition

Description: Essential for the control of the cell cycle at the G1/S (start) transition.

Antigen: Synthesized peptide derived from human Cyclin D1 around the phosphorylation site of Thr286

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: