

Anti-Human/Mouse/Rat Phospho-RPS6KA1 (Ser380) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01662

3th Edition

Description: Serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the transcription factor CREB.

Antigen: Synthesized peptide derived from human Rsk-1 around the phosphorylation site of Ser380

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: