

**Anti-Human/Mouse/Rat Phospho-PRKD2 (Ser876) Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA01669**

---

3th Edition

**Description:** Converts transient diacylglycerol (DAG) signals into prolonged physiological effects, downstream of PKC. Involved in resistance to oxidative stress.

**Antigen:** Synthesized peptide derived from human PKD2 around the phosphorylation site of Ser876

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Widely expressed.