

Anti-Human/Mouse/Rat Phospho-PRKD2 (Ser876) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01669

3th Edition

Description: Converts transient diacylglycerol (DAG) signals into prolonged physiological effects, downstream of PKC. Involved in resistance to oxidative stress.

Antigen: Synthesized peptide derived from human PKD2 around the phosphorylation site of Ser876

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed.