

Anti-Human/Mouse/Rat Phospho-CDC25A (Ser124) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01674

3th Edition

Description: Tyrosine protein phosphatase which functions as a dosage-dependent inducer of mitotic progression. Directly dephosphorylates CDK1 and stimulates its kinase activity. Also dephosphorylates CDK2 in complex with cyclin E, in vitro.

Antigen: Synthesized peptide derived from human Cdc25A around the phosphorylation site of S124.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: