

Anti-Human/Mouse Phospho-ERK 8 (Thr175/Y177) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01677

3th Edition

Description: In vitro, phosphorylates MBP.

Antigen: Synthesized peptide derived from human ERK 8 around the phosphorylation site of Thr175/Y177

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Widely expressed with a maximal expression in lung and kidney.