

Anti-Human/Mouse/Rat/Monkey Phospho-HNRNP K (Ser284) Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA01683

3th Edition

Description: One of the major pre-mRNA-binding proteins. Binds tenaciously to poly(C) sequences. Likely to play a role in the nuclear metabolism of hnRNAs, particularly for pre-mRNAs that contain cytidine-rich sequences. Can also bind poly(C) single-stranded DNA.

Antigen: Synthesized peptide derived from human hnRNP K around the phosphorylation site of S284.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: