

**Anti-Human CDC40 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02059**

---

3th Edition

**Description:** Associates with the spliceosome late in the splicing pathway and may function in the second step of pre-mRNA splicing.

**Antigen:** Recombinant human CDC40 protein

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**