

Anti-Human HMG20B Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02258

3th Edition

Description: Required for correct progression through G2 phase of the cell cycle and entry into mitosis. Required for RCOR1/CoREST mediated repression of neuronal specific gene promoters.

Antigen: Synthesized peptide derived from the N-terminal region of human BRAF35

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Ubiquitously expressed in adult tissues.