



**Anti-Human HMG20B Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02258**

---

3th Edition

**Description:** Required for correct progression through G2 phase of the cell cycle and entry into mitosis. Required for RCOR1/CoREST mediated repression of neuronal specific gene promoters.

**Antigen:** Synthesized peptide derived from the N-terminal region of human BRAF35

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Ubiquitously expressed in adult tissues.