

Anti-Human BRS3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02276

3th Edition

Description: The protein encoded by this gene is a G protein-coupled membrane receptor that binds bombesin-like peptides. This binding results in activation of a phosphatidylinositol-calcium second messenger system, with physiological effects including regulation of metabolic rate, glucose metabolism, and hypertension.

Antigen: Synthesized peptide derived from the Internal region of human BRS-3

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: In germ cells in testis. Lung carcinoma cells.