

Anti-Human/Mouse/Rat CABP7 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02379

3th Edition

Description:CABP7 (Calcium Binding Protein 7) is a Protein Coding gene. Among its related pathways are Visual Cycle in Retinal Rods. GO annotations related to this gene include calcium ion binding. An important paralog of this gene is CALN1.

Antigen:Synthesized peptide derived from the N-terminal region of human CaBP7

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: