

Anti-Human/Mouse/Monkey CCZ1/CCZ1B Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02724

3th Edition

Description:CCZ1 (CCZ1 Homolog, Vacuolar Protein Trafficking And Biogenesis Associated) is a Protein Coding gene. Among its related pathways are Vesicle-mediated transport and RAB GEFs exchange GTP for GDP on RABs. An important paralog of this gene is CCZ1B.

Antigen:Synthesized peptide derived from the C-terminal region of human CGI-43

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: