

Anti-Human/Mouse/Rat/Monkey CHEK2 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02762

3th Edition

Description:Regulates cell cycle checkpoints and apoptosis in response to DNA damage, particularly to DNA double-strand breaks. Inhibits CDC25C phosphatase by phosphorylation on 'Ser-216', preventing the entry into mitosis. May also play a role in meiosis. Regulates the TP53 tumor suppressor through phosphorylation at 'Thr-18' and 'Ser-20'.

Antigen:Synthesized peptide derived from human Chk2 around the non-phosphorylation site of T383.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:High expression is found in testis, spleen, colon and peripheral blood leukocytes. Low expression is found in other tissues.