

**Anti-Human/Mouse/Rat GJB5 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02779**

---

3th Edition

**Description:** This gene encodes a member of the beta-type (group I) connexin family. The encoded protein is a gap junction protein involved in intercellular communication related to epidermal differentiation and environmental sensing. This gene has been linked to non-small cell lung cancer.

**Antigen:** Synthetic peptide of human GJB5

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**