



Anti-Human Claudin 3 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02831

3th Edition

Description: Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.

Antigen: Synthesized peptide derived from human Claudin-3 around the non-phosphorylation site of Tyr219.

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: