

**Anti-Human /Mouse MAP2K3 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02833**

---

3th Edition

**Description:**Dual specificity kinase. Is activated by cytokines and environmental stress in vivo. Catalyzes the concomitant phosphorylation of a threonine and a tyrosine residue in the MAP kinase p38.

**Antigen:**Recombinant protein of human MAP2K3

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Abundant expression is seen in the skeletal muscle. It is also widely expressed in other tissues.