

Anti-Mouse/Rat/Monkey EIF4EBP1 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02890

3th Edition

Description:Regulates eIF4E activity by preventing its assembly into the eIF4F complex. Mediates the regulation of protein translation by hormones, growth factors and other stimuli that signal through the MAP kinase and mTORC1 pathways.

Antigen:Synthesized peptide derived from mouse 4E-BP1 around the non-phosphorylation site of T69.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: