

**Anti-Human/Mouse ADCY1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA02897**

---

3th Edition

**Description:** This gene encodes a form of adenylate cyclase expressed in brain. A similar protein in mouse is involved in pattern formation of the brain.

**Antigen:** Synthesized peptide derived from the Internal region of human A Cyclase I

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Brain, retina and adrenal medulla.