

Anti-Human /Mouse UNC5CL Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA02942

3th Edition

Description: This protein inhibits NF-kappa-B-dependent transcription by impairing NF-kappa-B binding to its targets. It interacts with p65/RELA and NFKB1. It has high levels expression in pancreas, liver and kidney. It contains one death domain and one ZU5 domain.

Antigen: Synthetic peptide of human UNC5CL

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in pancreas, liver and kidney.