

Anti-Human DDIAS Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA06497

3th Edition

Description:DDIAS (DNA Damage Induced Apoptosis Suppressor) is a Protein Coding gene.

Antigen:Synthesized peptide derived from the Internal region of human NOXIN.

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Highly expressed in colorectal and lung cancer tissues.