

**Anti-Human DDX53 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA09270**

---

3th Edition

**Description:** This intronless gene encodes a protein which contains several domains found in members of the DEAD-box helicase protein family. Other members of this protein family participate in ATP-dependent RNA unwinding.

**Antigen:** Recombinant protein of human DDX53

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in testis. Wide expression in various cancer tissues and cancer cell lines.