

Anti-Human CEACAM21 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11728

3th Edition

Description:CEACAM21 is a gene which encodes carcinoembryonic antigen.It belongs to the immunoglobulin superfamily-CEA family. Three isoforms of the human protein are produced by alternative splicing.High levels of expression in hematopoietic system, brain and muscle.

Antigen:Recombinant protein of human CEACAM21

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: