

**Anti-Human CEACAM21 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA11728**

---

3th Edition

**Description:**CEACAM21 is a gene which encodes carcinoembryonic antigen.It belongs to the immunoglobulin superfamily-CEA family. Three isoforms of the human protein are produced by alternative splicing.High levels of expression in hematopoietic system, brain and muscle.

**Antigen:**Recombinant protein of human CEACAM21

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**