

Anti-Human/Mouse RNF7 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA11764

3th Edition

Description:The protein encoded by this gene is a highly conserved ring finger protein. It is an essential subunit of SKP1-cullin/CDC53-F box protein ubiquitin ligases, which are a part of the protein degradation machinery important for cell cycle progression and signal transduction. This protein interacts with, and is a substrate of, casein kinase II (CSNK2A1/CKII). The phosphorylation of this protein by CSNK2A1 has been shown to promote the degradation of I κ B α (CHUK/IKK- α /IKBKA) and p27Kip1(CDKN1B).

Antigen:Recombinant protein of human RNF7

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1 μ g/ml.

Specificity:Expressed in heart, liver, skeletal muscle and pancreas. At very low levels expressed in brain, placenta and lung.