

**Anti-Human/Mouse FMNL1 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12083**

---

3th Edition

**Description:** This gene encodes a formin-related protein. Formin-related proteins have been implicated in morphogenesis, cytokinesis, and cell polarity. An alternative splice variant has been described but its full length sequence has not been determined.

**Antigen:** Recombinant protein of human FMNL1

**Form:**

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas.