

**Anti-Human/Mouse/Rat FKBP8 Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12086**

---

3th Edition

**Description:**The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. Unlike the other members of the family, this encoded protein does not seem to have PPIase/rotamase activity. It may have a role in neurons associated with memory function.

**Antigen:**Recombinant protein of human FKBP8

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Widely expressed. Highest levels seen in the brain.