

Anti-Human/Mouse/Rat ENTPD5 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12143

3th Edition

Description:The protein encoded by this gene is similar to E-type nucleotidases (NTPases)/ecto-ATPase/apyrases. NTPases, such as CD39, mediate catabolism of extracellular nucleotides. ENTPD5 contains 4 apyrase-conserved regions which is characteristic of NTPases.

Antigen:Recombinant protein of human ENTPD5

Form:

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Expressed in adult liver, kidney, prostate, testis and colon. Much weaker expression in other tissues.