

**Anti-Human/Mouse/Rat FSHR Polyclonal Antibody**

**Polyclonal Antibody**

**Cat.NO.: PA12206**

---

3th Edition

**Description:**The protein encoded by this gene belongs to family 1 of G-protein coupled receptors. It is the receptor for follicle stimulating hormone and functions in gonad development. Mutations in this gene cause ovarian dysgenesis type 1, and also ovarian hyperstimulation syndrome. Alternative splicing results in multiple transcript variants.

**Antigen:**Recombinant protein of human FSHR

**Form:**

**How to use:**1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:**PBS (pH7.4) containing 1% BSA

**Application:**This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:**Sertoli cells and ovarian granulosa cells.