

Anti-Human DDX53 Polyclonal Antibody

Polyclonal Antibody

Cat.NO.: PA12231

3th Edition

Description: This intronless gene encodes a protein which contains several domains found in members of the DEAD-box helicase protein family. Other members of this protein family participate in ATP-dependent RNA unwinding.

Antigen: Recombinant protein of human DDX53

Form:

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Expressed in testis. Wide expression in various cancer tissues and cancer cell lines.